DIALTANC® / CITRATANC

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The product goes directly from the tank of the circuit to the dialysis machine.

There is no handling from the user.

Better stability throughout the dialysis session.

Alarm and informational panel from the circuit in the dialysis room.

Bacterial contamination safety system.

The tanks are connected to a pumping system which supplies the concentrate directly into the dialysis machines.

The product is to be used in centres which have an online distribution system with sanitary compatible materials.

ADVANTAGES:

- It reduces product waste.
- It is an easy-to-use, safe and aseptic system.
- It saves time in preparation.
- It avoids moving canisters around the room.
- It reduces plastic and cardboard waste.

TECHNICAL CHARACTERISTICS:

DIALTANC® is available in acid concentrate (with acetic acid 1/35, 1/36, and 1/45 dilution rate) and bicarbonate.

CITRATANC is available in acid concentrate (with citric acid, 1/35, 1/36, and 1/45 dilution rate).

The acid concentrate is packed directly into the tank.

The bicarbonate concentrate is introduced into a one use blown polyethylene bag placed inside the tank. The concentrate goes directly from the bag to the distribution circuit to the machine. As the product is used, the bag folds back in itself. By using a vacuumed bag the contact of the product with the exterior is avoided, any possible contamination is prevented and the characteristic of the product, pH and conductivity, are preserved.

The product is always sealed.

DIALTANC® / CITRATANC have a wide range of formulations, allowing choosing the best treatment for each patient.

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TANK'S CHARACTERISTICS:

The tank has a capacity of 1000 or 600 litres, it is clear, made of high density polyethylene stabilized with ultraviolet (UV) rays. The body is placed inside a galvanized tube support and fitted into a plastic or wood pallet (stackable).

The dimensions are: for 1000 litres tank, 1000 mm x 1200 mm x 1170 mm, weight, 59 to 71kg; for 600 litres, 800 mm x 1200 mm x 995 mm, weight, 46 to 58kg.

The tank has 2 outlets, one on the top part of it which has to be loose when starting to use and the other one, on the lower part, where the circuit is connected for its distribution.